

DynamiQ display

For immediate display of analysis results and instrument status


Unlock the potential of the optional, external, explosion-proof DynamiQ display. Experience the convenience of immediate viewing of results, instrument status and direct GC control. Be empowered to elevate gas analysis performance by using the DynamiQ display.



Connection

The DynamiQ display can be connected to any Ex e-electrical junction box with the communication cable, which can be placed in the same explosive hazardous zones as the DynamiQ-X micro-GC.

This enables users to use their own cables, e.g. an ethernet cable and separate power cable that connects to the GC.



Facts

Your benefits

- Easy and secured viewing of GC results and instrument status
- Direct GC control

Use cases

- Immediate analysis of GC results
- Immediate access to GC control in Ex operation areas

Easy to use

- Easy navigation through the different screens of the display
- No prior knowledge of the DynamiQ software required